

Following are simple, step-by-step instructions on how to set-up the Dolphin Wireless Timing System for a Meet. For more detailed instructions refer to the information in the [CVS Dolphin User Guide](#).

### **Part One – Set-up the timing equipment**

1. Turn on the **Stopwatches** (Press and hold the RESET button).
2. Mount the **Dolphin Infinity** starting device on the tripod.
3. Connect the **Microphone** to the Dolphin Infinity starting device.
4. Connect the **Dolphin Start System** to the Dolphin Infinity starting device.
5. Connect the **Dolphin Transceiver** to the CVS-DOLPHIN computer.
6. Turn on the **Dolphin Start System** (Press and hold the RESET button).
7. Turn on the **Dolphin Infinity** starting device.

### **Part Two – Set-up the computers**

1. Connect the **CVS-DOLPHIN** computer to the **CVS-MEETMANAGER** computer using the yellow ethernet cable.
2. Connect the printer to the **CVS-MEETMANAGER** computer.
3. Turn on the printer.
4. Turn on the **CVS-DOLPHIN** computer. Login with credentials - CVSwimming/CVSwimming
5. Turn on the **CVS-MEETMANAGER** computer. Login with credentials – CVS/CVS
6. **CVS-MEETMANAGER:** In File Manager, select **CTSDolphin (\\CVS-DOLPHIN) (M:)**
7. **CVS-MEETMANAGER:** Delete all \*.D03 and \*.D04 files on CTSDolphin (\\CVS-DOLPHIN) (M:)
8. **CVS-MEETMANAGER:** Start the **Meet Manager** program and load the Swim Meet database.
9. **CVS-DOLPHIN:** Start the **Dolphin.exe** program.

### **Part Three – Perform a Test Start and get times into the Meet Manager database**

1. Perform a Test Start. After 30 seconds, press **STOP** on all stopwatches.
  2. When all stopwatches have been stopped, the CTS Dolphin window will show a placing for each lane.
  3. **CVS-DOLPHIN:** Press **RESET TIMERS("r")** on CTS Dolphin window. This will commit the times to the Dolphin program, reset the stopwatches and ready the Dolphin Start System for the next race..
  4. **CVS-MEETMANAGER:** Meet Manager → Run → Interfaces → Timer (CTSS) Pool 1 → Select Data Set stored from CTSS.
  5. **CVS-MEETMANAGER:** Press **Update Data Location**
  6. **CVS-MEETMANAGER:** Select CTSDolphin (\\CVS-DOLPHIN) (M:) → **Open**.
  7. **CVS-MEETMANAGER:** Press **Next Meet** until a Meet with the time and date of the Test Start appears → **Close**.
  8. **CVS-MEETMANAGER:** Press **Get Times : F3**
  9. **CVS-MEETMANAGER:** Select file \*-0001.d04 → **OK**
-